

# Prepping libraries for sequencing submission

```
source("../genomics/scripts/lab_helpers.R")

##
## Attaching package: 'dplyr'
## The following objects are masked from 'package:stats':
##
##   filter, lag
## The following objects are masked from 'package:base':
##
##   intersect, setdiff, setequal, union
lab <- write_db("Laboratory")

## Loading required package: DBI
print(params)

## $date
## [1] "2018-09-13"
##
## $facility
## [1] "PrincetonLSI"
##
## $first_baits
## [1] "B001"
##
## $frag_length
## [1] 430
##
## $last_baits
## [1] "B003"
##
## $library_type
## [1] "ddRADSeq"
##
## $num_lanes
## [1] 2
##
## $pe_or_se
## [1] "se"
##
## $platform
## [1] "hiseq"
##
## $seq
## [1] "SEQ18"
# pull in the baits that are going to be pooled for sequencing
baits <- lab %>%
  tbl("baits") %>%
  collect()
```

```

# update the db to reflect pooling
change <- baits %>%
  filter(baits_id >= params$first_baits & baits_id <= params$last_baits) %>%
  mutate(seq = params$seq)

baits <- change_rows(baits, change, "baits_id")
# dbWriteTable(write_lab, "baits", dat_db, row.names=F, overwrite = T)

# calculate the percent concentration of DNA in each bait
change <- change %>%
  mutate(vol = 20,
    perc = as.numeric(ng_ul)/as.numeric(vol))

# calculate the percent concentration of the pooled baits
tot <- change %>%
  # sum of percents
  summarise(tot_perc = sum(perc)) %>%
  # sum of volumes
  mutate(vol = change$vol[1] * 3,
    # calculate the ng_ul of the combined DNA
    ng_ul = tot_perc * vol,
    # calculate the nM of the combined DNA (430bp fragment length)
    nM = (ng_ul/(660*params$frag_length))*10-6,
    # calculate the volume of library to produce 10nM
    seq_vol = (30*10)/nM) # 30 because that is our desired final volume to send to seq, 10 because we w

knitr::kable(tot) %>%
  kableExtra::kable_styling()

```

```

## Warning in kableExtra::kable_styling(.): Please specify format in kable.
## kableExtra can customize either HTML or LaTeX outputs. See https://
## haozhu233.github.io/kableExtra/ for details.

```

tot_perc	vol	ng_ul	nM	seq_vol
0.1074	60	6.444	22.70613	13.21229

```

# add to seq table in the db
dat_db <- lab %>%
  tbl("sequencing") %>%
  collect()

# create a row to add to the seq table in the db
change <- tibble(seq_id = params$seq, date = params$date, library_type = params$library_type, facility =

# -----
# add columns to seq db (only have to do this once today)
dat_db <- dat_db %>%
  mutate(vol_pooled_libs = NA, vol_TE = NA)
dat_db <- dat_db %>%
  select(seq_id:contents, vol_pooled_libs, vol_TE, library_type:notes)
# -----

```

```

dat_db <- rbind(dat_db, change) %>%
  mutate(vol_pooled_libs = formatC(vol_pooled_libs, width = 1),
         vol_TE = formatC(vol_TE, width = 1))

# dbWriteTable(lab, "sequencing", dat_db, row.names=F, overwrite = T)

```

Paste the qubit results into the PCR page of the Sample\_Data file as directed above in the qubit protocol.

On the PCR page, in the  $\mu\text{L}$  PCR product to make  $30\mu\text{L}$   $10\text{nM}$  column, pipet that number into the sequencing vial. Then pipet the number from the water column.

Label a vial with the SEQ number on top (SEQ 04) and SEQ # date nM ng/ $\mu\text{L}$  indices on the side for example (this really does fit):

Pinsky SEQ 04 12-10-2014 10 nM Indices 4,5,6,7

Stuff a kim wipe into a  $50\text{mL}$  falcon tube, place the diluted  $10\text{nM}$  sequencing tube inside, stuff in another kim wipe, screw on the cap and place in a styrofoam cooler with an icepack. We want to keep the sample cool but we do not want to freeze it. If it isn't ready to depart immediately, keep the  $50\text{mL}$  falcon with the sequencing tube in it in the fridge.

Place the original pooled sequencing sample vial in the freezer box in the  $-20$

## Update 4/13/2015 - MRS

It is possible to calculate the concentration of each pool and then add  $100\text{ng}$  of each:

Based on the qubit quantification (X), determine the volume of sample (Y) that will supply  $100\text{ng}$  to the sequencing pool:

$$100\text{ng} \div X \text{ ng}/\mu\text{L} = Y \mu\text{L}$$

The concentration of this pooled sequencing sample is calculated in the next cell by dividing  $400$  by the sum volume ( $400$  because you added  $100\text{ng}$  of  $4$  pools) - this value should calculate automatically in the spreadsheet

We are using  $430\text{bp}$  fragments ( $375$  from pippin plus  $55$  added by illumina primer for index and flowcell attachment) - enter this value into the spreadsheet

The spreadsheet will calculate  $\text{nmol}/\mu\text{L}$  by dividing the concentration by the product of  $660$  and the number of basepairs ( $430$ ) ...  $660 * 430 = 283,800$

The spreadsheet will calculate  $\text{nM}$  (aka  $\text{nmol}/\text{L}$ ) by multiplying  $\text{nmol}/\mu\text{L}$  by one million

Calculate the  $4$  pools independently on the PCR sheet in the Sample\_Data file

Sum the volume on the Sequencing sheet in the Sample\_Data file

Combine the  $4$  pools based on the volume calculated

The next two cells in the Sequencing spreadsheet calculate the volume of sample and the volume of water to be combined in this new tube - the volume of pooled sequencing sample is calculated by dividing  $300$  by the  $\text{nM}$  and the volume of water is calculated by subtracting the volume of sequencing sample from  $30$ .

Combine the sequencing sample and Pure water in the new tube.

Stuff a kim wipe into a  $50\text{mL}$  falcon tube, place the diluted  $10\text{nM}$  sequencing tube inside, stuff in another kim wipe, screw on the cap and place in a styrofoam cooler with an icepack. We want to keep the sample cool but we do not want to freeze it. If it isn't ready to depart immediately, keep the  $50\text{mL}$  falcon with the sequencing tube in it in the fridge.

Place the original pooled sequencing sample vial in the freezer box in the -20

If sending to Princeton, don't generate the PO until the seq is complete (particularly if you are submitting more than one seq within a month). Ask Princeton if they can add some description to the quote to make it look different than the rest of the quotes they have sent (like add the seq number). We have had problems in the past with multiple quotes from similar dates with no other difference on the page. Had to add a line of text saying "This invoice is for sequencing done to APCL12 DNA on February 6th, 2016".

If sending the sample to Princeton, Malin will fill out the submission form on this webpage: <https://htseq.princeton.edu>

Princeton also has guidelines on this webpage: <http://www.princeton.edu/genomics/sequencing/instructions/sample-preparation/>

Their guidelines differ from our personal contact's recommendations at Princeton and we will follow our own set of guidelines unless instructed otherwise by Malin.